

# VELOX<sup>TM</sup>

## Solvent Waste System

### Velox Heavy Duty Solvent Waste System

Zodiac Life Sciences introduces Velox Heavy Duty HPLC Solvent Waste system. This includes HD carboys, exhaust filters, cap kits. These provide efficient and affordable solutions for waste disposal. Each system is manufactured to safely collect waste from one or more HPLC systems while the carbon exhaust filter improves air quality by capturing potentially harmful vapors emitted from solvents.



Carbon exhaust filter traps solvent vapors



90° Adaptor included for maximum storage flexibility



Tube fittings secure waste tubes to carboy providing a closed system that prevents vapors from escaping

### Ordering Information

Item Description	Part Number
Velox Solvent Waste System (10 Ltrs)	ZLS.332-3413
Velox Solvent Waste System (20 Ltrs)	ZLS.332-4413