

Velox Solvent Safety Cap



Advantages:

The closed safety system of Zodiac offers a calming protection and ensures an economical work in laboratory. The new, improved air valve combines valve and filter functions. As usual, ventilation occurs during removal; harmful solvent vapors are blocked. At the same time, the valve membrane absorbs dust and contaminant particles from the incoming air. The valve also fits your existing Velox Safety Caps without any technical modification.



No more harmful vapors

Thanks to the integrated air valve.

- No eluent contamination
 Receptacles remain securely closed so the results of your analysis stay correct.
- No tube slippage
 No accidental intake of air into the HPLC system.
 - No interruption of analytical and work processes due to air pockets in the tube.
- Easy container changes
 Thanks to a freely-rotatable cap, even with tubes installed. Without twisting or "tangled tubes".
- Confidence during audits
 With Velox Solvent Safety Cap you pass quality and safety inspections with confidence.

Air Valve:

- The new, improved air valve combines valve and filter functions. As usual, ventilation occurs during removal.
- Harmful solvent vapours are blocked. At the same time, the valve membrane absorbs dust and contaminant particles from the incoming air.
- The valve also fits your existing Velox Safety Caps without any technical modification.





Sign-in the installation date:

Never miss a chance again

The air valve with the part number is delivered with a change-advice-label.

Just sign-in the current date at the installation and you are safe in observing the service life.

Ordering Information:

Part No.	ZLSIVSC01	ZLSIVSC02	ZLSIVSC03	ZLSIVSC04	ZLSIVSC05	ZLSVSCS01	ZLSAV
Description	Safety Cap I	Safety Cap II	Safety Cap III	Safety Cap IV	Safety Cap VI	Safety Cap Set	Air Valve
Thread size	GL-45	GL-45	GL-45	GL-45	GL-45	GL-45	GL-45
Connector 3.2mm Ø OD	1X	2X	3X	4X	6X	4X	-