






## Velox Solvent waste systems

## Introduction:

Fluids in HPLC waste containers are overwhelmingly dangerous containing potent volatile organic compounds (VOC) posing potent risk to health of Lab workers & lab safety.

Zodiac came with an innovative solution of introducing **VELOX Solvent waste system** with activated carbon filters not only creating an organised containment for HPLC solvent waste but also providing an efficient way for collection and disposal of volatile vapours.



## Advantages

Filters system against harmful vapours


- Capillary fixations with PFA Fittings
- Several systems can be connected to 1 container
- Different size containers compatible to any space condition



90 Adaptor included for maximum storage flexibility



Normal cap is additionally provided for convenient transport of the wastes





Blind fittings are provided to block the remaining ports when not in use

## ORDERING INFORMATION

Item Description	Part Number
Velox Solvent Waste System - 5Ltrs	ZLS.332.2419
Velox Solvent Waste System - 10Ltrs	ZLS.332-3413
Velox Solvent Waste System - 20Ltrs	ZLS.332-4413